# CHEM 450 Expt. 2: DETERMINATION of CALCIUM in "URINE" SAMPLES by VISIBLE (Vis) SPECTROSCOPY

Advanced Preparation: To be completed before coming to the lab.

Show calculations to your instructor before preparing any of the solutions below.
100 mL of 0.00100 M Ca<sup>2+</sup> from 400 ppm Ca<sup>2+</sup> stock solution
200 mL of a 1.00 x 10<sup>-3</sup> M aqueous solution of Arsenazo III (FW 820.33 g/mole)
1 L of 0.1 M Imidazole buffer from the hydrochloride salt, adjusted to pH 6.5 with either 0.1 M HCl or 0.1 M NaOH

## Background

Soluble calcium (as  $Ca^{2+}$  ions) in aqueous samples is usually determined by EDTA titration using Eriochrome Black T as indicator. In very dilute solutions, titration becomes impractical to use due to low levels of calcium ions which make the color change at the equivalence point difficult to see.

Calcium is usually measured to screen for or monitor diseases of the bone or calcium regulation disorders (that is, diseases of the parathyroid gland or kidneys). Urinary calcium levels aid the clinician in understanding how the kidneys handle calcium in certain disease states involving calcium balance including diseases of the parathyroid gland. Urinary calcium levels are also essential in the medical evaluation of kidney stones. (http://medlineplus.gov/) Urine calcium is usually measured in a sample taken from all the urine produced in a 24-hour period. Normal results may vary from lab to lab. Test results are affected by the amount of calcium in the diet. Therefore, your health professional will consider the amount of calcium in your diet when interpreting your urine calcium levels. (http://www.webmd.com/hw/lab\_tests/hw27965.asp.)

Calcium in urine				
Low amount of calcium in diet:	Less than 150 milligrams (mg)/24-hour sample			
Average amount of calcium in diet:	100–250 mg/24-hour sample			
High amount of calcium in diet:	250–300 mg/24-hour sample			

Ultraviolet/visible (UV/vis) absorption spectroscopy is a technique that can be used to determine low levels of calcium. Because aqueous solutions of calcium salts are colorless and therefore do not absorb in the near UV and visible region, an organic complexing agent will be used in this determination. Arsenazo III (1,8-dihydroxynaphthalene-3,6-disulphonic acid-2,7-bis[(azo-2)phenylarsonic acid]) forms highly colored 1:1 complexes with calcium and other metallic ions, allowing calcium levels to be determined using UV/vis spectroscopy. Calcium-arsenazo complex is blue or purple in color, depending on pH, while uncomplexed arsenazo is wine-red. This complexation process is also utilized in the determination of calcium ions in serum and urine. Calcium measurement in biological samples is important in the treatment of a variety of bone diseases and chronic renal disease.



Because the complexation between  $Ca^{2+}$  ions and arsenazo is pH dependent, the solutions and samples will be buffered at pH 6.5.

## Chemicals Used

 $1.00 \ge 10^{-3}$  M Arsenazo III, disodium salt  $1.00 \ge 10^{-3}$  M Ca<sup>2+</sup> from 400. ppm Ca<sup>2+</sup> stock pH 6.5 Imidazole buffer (see advanced prep) One 24-hr "urine" sample

## Procedure

## A. Preparation of Arsenazo and Calcium Chloride Solutions

1. (a) Using the FW of Arsenazo III (820.33 g/mole for the disodium salt) calculate the amount in grams needed to prepare 200.0 mL of a 1.00 x  $10^{-3}$  M aqueous Arsenazo III solution. Enter your calculation in your lab notebook. *You do not have to weigh this. Instead, proceed to step 1(b).* 

(b) Calculate the volume of 400. ppm  $Ca^{2+}$  stock solution needed to prepare 50.00 mL of 1.00 x  $10^{-3}$  M  $Ca^{2+}$ . Calculate the equivalent ppm concentration of this new dilution. Record <u>both</u> concentrations in your notebook.

2. *For each group*: Prepare 50.00 mL of the  $1.00 \times 10^{-3} \text{ M Ca}^{2+}$  solution following your calculation above. The Arsenazo solution is already prepared based on the amount calculated above and is ready for use.

# B. Preparation of 0.1 M pH 6.5 Imidazole Buffer

1. Calibrate a pH meter using pH 7.00 and 10.00 buffers.

2. Calculate the amount of imidazole hydrochloride salt needed to prepare 100.0 mL of 0.10 M solution. Weigh the calculated amount on a weigh boat and transfer quantitatively (*ask your instructor if you don't know what this means*) into a 200-mL beaker. Dissolve in ~80mL of water. Check the pH of your buffer. Is it below or above the desired pH?

3. With constant stirring, add either 3.0 M NaOH or 3.0 M HCl dropwise until pH 6.5 is reached. Transfer into a 100-mL volumetric flask. Dilute to the mark with DI water. Cap and shake well.

# C. Preparation of Calcium-Arsenazo III Standard Solutions and Unknown.

1. Obtain eight 25-mL volumetric flasks and label according to Table 1.1 Flask ID.

2. Obtain one unknown "urine" sample from the hood for calcium determination. Write down your unknown number.

3. Quantitatively transfer 5.00 mL of your unknown into a 50.00 mL volumetric flask. Dilute to the mark with deionized water. (*Save the remaining unknown sample for a future lab*)

4. Label a 25-mL volumetric flask "Patient (Unknown #)" and another 25-mL volumetric flask 25-mL volumetric flask "Patient (Unknown #) Duplicate".

5. Transfer 1.00 mL unknown solution from step C3 into each of the labeled 25.00-mL volumetric flasks using a calibrated micropipet. *Do not dilute to the mark yet.* Proceed to the next step.

6. Following Table 1.1 prepare a <u>solvent blank solution</u>, <u>four dilutions</u> of the Ca-Arsenazo standard solution from part A buffered to pH 6.5 (called <u>working standards</u>), a solution of <u>Arsenazo III</u>, and your <u>unknown urine samples</u> by combining the specified reagents into a 25.00-mL volumetric flask.

7. Dilute each flask to the 25.00-mL mark with deionized (DI) water. Cap and shake well. Allow at least 15 minutes for the complexation to take place before taking any absorbance measurements.

8. Using both the ppm and moles/L concentrations of your stock Ca-Arsenazo standard solution from part A, calculate the adjusted concentration of  $Ca^{2+}$  in the Ca-Arsenazo *working standards* following dilution with DI water to 25.00 mL. Enter these concentrations in the last two columns of Table 1.1. [*Hint:* In all dilutions  $Ca^{2+}$  is the limiting reagent and Arsenazo III forms a 1:1 complex with  $Ca^{2+}$ ].

Table 1.1. Amounts of reagents needed to	prepare Ca <sup>2+</sup> - arsenazo III standard solutions.
0	1 1

Flask ID	Volume of	Volume of	Volume	Total	Concentration	Concentration of
	$1.00 \ge 10^{-3} M$	$1.00 \times 10^{-3} M$	of	Volume	of Ca <sup>2+</sup> in the	Ca <sup>2+</sup> in the
	Arsenazo III	Ca <sup>2+</sup>	buffer	(mL)	working	working standard
	(mL)	(mL)	(mL)		standard soln.	soln.
					(mol/L)	(ppm or mg/L)
Blank	0	0	10.00	25.00		0
Arsenazo III	1.00	0	10.00	25.00		0
Std 1	1.00	0.250	10.00	25.00		
Std 2	1.00	0.500	10.00	25.00		
Std 3	1.00	0.750	10.00	25.00		
Std 4	1.00	1.00	10.00	25.00		
Patient	1.00	0	10.00	25.00		
Patient	1.00	0	10.00	25.00		
duplicate						

# D. Spectroscopic Determination of the Amount of Ca<sup>2+</sup> in the Unknown

Follow the operating instructions for the double-beam UV/Vis spectrophotometer.

1. Obtain the absorption spectrum of the blank solution.

2. Using the visible range of 450 to 750 nm wavelength obtain the absorption spectrum of your Arsenazo solution. Compare this spectrum with that of the Ca-Arsenazo complex by overlapping the spectrum of Std 3 with that of Arsenazo. Print and label each spectrum accordingly.

*Question:* From the result of the overlapped spectra, what wavelength should be used to quantify  $Ca^{2+}$  in the Ca-Arsenazo III complex and the unknowns? Why?

3. Depending on the display capability of the spectrophotometer, you may have to clear the absorption spectra obtained above before proceeding with the next steps.

4. Following the instructions on "Quantification" (also Quantitation) and using the wavelength of maximum absorption,  $\lambda_{max}$ , generate a Calibration Curve by acquiring the spectra of the four working standards using an expanded wavelength range of 600 to 700 nm and plotting the resulting Absorbance vs. Concentration data. Select the option that displays the line equation and correlation coefficient.

5. Measure the absorbance of your unknown solution(s) and, using your calibration curve, calculate the concentration of  $Ca^{2+}$  in your unknown. Remember to use the dilution factor (all your dilutions multiplied together) when calculating the concentration of calcium ion in your original, undiluted unknown solution(s).

6. Print your data for inclusion in your laboratory report.

#### Treatment and Reporting of Data

1. Calculate the mean ppm  $Ca^{2+}$ , standard deviation, and % relative standard deviation (RSD) of *duplicate* measurements for your unknown. Report these in your RESULTS with the unknown number.

2. Is the  $[Ca^{2+}]$  in your unknown normal or high? Discuss in enough details the implications of abnormal (both high and low) levels of calcium in urine. Cite your references.

#### Reference

"Photometric Determination of Micro Amounts of Calcium with Arsenazo III." *Anal. Chim. Acta*, 53 (1971) 194-198.